

# Protekt Akra Bariatric™ Bed Series



The **Protekt Akra Bariatric™ Bed Series** offers the lightest 42" bed on the market. It's extremely easy to maneuver, transport, and assemble. Weighing in at only 200 lbs. takes nothing away from the Akra's independently tested strength and durability. Available in 42", 48", and 54" widths provide patients with extra space and comfort to support their recovery.

## FEATURES

- Ultra lightweight frame design offers providers easier and more manageable bed deliveries and assemblies
- 42" bed ships complete with rails in only 2 boxes
- Ultra lightweight frame is easily set up by one individual
- Full electric bariatric homecare bed with dual motor assembly for head and frame adjustability
- Square tube steel frame slat deck provides durable strength, no-snag mattress support and is easy to clean
- Scratch resistant wood grain panel bed ends include hi/lo motors for frame height adjustability
- UL approved motor is self-contained for reduced weight and noise and can be easily installed & removed with patient in bed
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Rail accessories meet FDA and ASTM entrapment guidelines and standards

## SPECIFICATIONS

	PMBAB42 - 42"	PMBAB48 - 48"	PMBAB54 - 54"
<b>Deck Height Range:</b>	13" min.	13" min.	13" min.
	26" max.	24.5" max.	24.5" max.
	* Removing casters will lower range 6"		
<b>Sleep Surface:</b>	80"L x 42"W	80"L x 48"W	80"L x 54"W
<b>Overall Bed:</b>	87.5"L x 42"W	88.6"L x 48"W	88.6"L x 54"W
<b>Weight Capacity:</b>	600 lbs.	750 lbs.	1000 lbs.
<b>Product Weight:</b>	200 lbs.	298 lbs.	307 lbs.
<b>Warranty:</b>	Limited Lifetime on Frame & Welds		
	5 Years on Motor & Hand Controls		
	1 Year on All other Parts & Components		

Item #	Description	HCPCS Code
PMBAB42	Bariatric Full-Electric Bed w/Half Rails 42"	EO301
PMBAB48	Bariatric Full-Electric Bed w/Half Rails 48"	EO302
PMBAB54	Bariatric Full-Electric Bed w/Half Rails 54"	EO302
BP-HRBA	Half Length Side Rails f/Bariatric Bed (for 42"/48"/54")	EO305

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com