


Protekt Akra-FE™ Full-Electric Bed

The **Akra-FE™** is multi-functional, height adjustable and designed for user and caregiver convenience. In the event of a power outage, a 9V battery backup can be used to allow movement to the lowest position. The sturdy lightweight frame and snap on or off low-voltage DC motor box ensure quick, efficient set-ups and easy routine maintenance.



Half Rails



Full Rails



Pendant

FEATURES

- The entire bed, including the rails comes complete in only (2) boxes
- Quick & simple motor attachment or removal for service or cleaning
- Welded steel split frame construction
- Easy-to-use control adjusts upper body, knees, and height of bed
- Head and foot sections can be raised independently or simultaneously
- Large button, simple to use hand pendant
- Scratch-resistant bed ends keep the bed looking new
- 9V battery back-up (battery not included)
- Headspring and foot assembly are recessed to help reduce pinching and injuries
- Beds include 2 locking and non-locking casters
- Quiet, efficient ETL-approved DC low voltage motors

SPECIFICATIONS

Deck Height:	15" min. no castors 23" max. w/castors
Sleep Surface:	80"L x 36"W
Overall Bed:	86"L x 36"W
Weight Capacity:	450 lbs.
Warranty:	Limited Lifetime Frame 2 Years on Motor

SHIPPING DIMENSIONS

Component	Carton Size	Shipping Weight
Foot/Head Springs, Motors, Pendant	47"x35.5"x9"	112 lbs. w/HR 117 lbs. w/FR
Head/Foot Bed Ends, Casters	38"x26"x5"	57 lbs.

Item #	Description	HCPCS Code
PMFEBHR	Full-Electric Bed w/Half Rails (No Support Surface)	E0266
PMFEBFR	Full-Electric Bed w/Full Rails (No Support Surface)	E0266
PMFEBHRM	Full-Electric Bed w/Half Rails w/Protekt 100 Foam Mattress (Group 1)	E0265
PMFEBFRM	Full-Electric Bed w/Full Rails w/Protekt 100 Foam Mattress (Group 1)	E0265
PMFEB	Full-Electric Bed w/OUT Rails (No Support Surface)	E0297
BP-HR	Half Length Side Rails (pr)	E0305
BP-FR	Full Length Side Rails (pr)	E0310